# dental diamant *dental diamond*

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## Diamonds made by Fis - the alternatives

*We have been manufacturing diamond tools for dental as well as diamond and CBN tools for industrial applications since more than 75 years. Shapes, sizes and abrasive layers of our products are very carefully matched for optimum use in all applications. Custom products have been manufactured for many of our customers. This has resulted in "smart" solutions which are unparalleled in solving technological requirements, offering high quality products.*

### »High-tech is a standard«

*For dental applications, our production ranges from FG burs to plaster trimmer wheels. We manufacture our FG diamonds in varying grit sizes between 504 and 554 (acc. to ISO standards), with choice of one of two types of grit: natural and/or synthetic (S) diamond (limited to certain sizes and types). The natural diamond crystals typically have an irregular shape and sharp cutting edges which ensure a quick material removal. Layers containing the blocktype, synthetic diamond crystals provide a professional preparation with less surface roughness. Those FG diamonds have a black coating to increase the optical contrast between layer and tooth. Only hardened, stainless chromium/nickel steels are used as basic material for the "fis" diamond tools (except a few cutting discs). The tolerances of the shaft diameters not only meet the requirements of the ISO standard, but exceed them. With quality in mind, we grind the shafts to close tolerances for precise fit in the handpiece chuck which minimizes runout error and protects the chuck. In addition, we have introduced a new generation of diamonds (HighQPlus). These multilayer FG tools provide long life in rough cutting applications.*

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## Directions for use of rotary diamond instruments (FG and right angle type)

The following instructions provide best results in preparations through careful handling of the diamond instruments and their application. Consideration of the teeth as well as perfect hygienic conditions are the main concern.

- In terms of cleaning/sterilization please follow the information and instructions for cleaning and sterilizing electroplated dental diamond burs with FG or right angle shank on page 4.
- Instruments are to be selected by the head diameter to prevent from damaging neighbour teeth or tissue.
- Rounded shapes should be favoured.
- Instruments with coarse grit are not always indicated; usually they require a finish grinding operation.
- Trays or original packing should be used for storing.
- The shafts have to be secured properly and checked for minimum runout before use.
- The burs should be set to the working speed according to the graphical diagram. It should be set before applying the bur to the object.
- Spray cooling with a minimum flow of 50 ml/min must be provided. Additional cooling is required as burs longer than 21 mm or heads bigger than 2 mm are being used.
- Forces should not exceed 2 N (appr. 200 p). Higher forces do not result in a quicker material removal, but increase the thermal load on the tooth.
- Worn or damaged burs must be replaced immediately.
- Burs should be cleaned in a "cleaning solution" first after use. See manufacturer's instructions for details.
- Burs can be brushed to remove material sticking to the diamond layer.
- Dental diamond instruments are made of stainless materials. However, they must not get in contact with acids.
- Our dental diamond tools are manufactured from high-quality components providing a long life time. However we recommend not to use the instruments after an expiration time of 10 years.

## Tools with sintered bond

Our tools with sintered bond are particularly suitable for processing ceramic materials and metal bondings. The high degree dimensional stability and the long service life are ensured by a special production technique.

## Directions for use of rotary diamond instruments (Lab tools with sintered bond)

Tools with diamond grit in a sintered bond offer a variety of applications combined with a long tool life. Some simple rules ensure maximum use of their advantages.

- The tools should be set to the working speed according to the graphical diagram.
- A slightly oscillating motion provides low cutting forces by a short length of contact.

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
- Forces should be kept low to avoid heating and excessive tool wear.
- Cutting discs with sintered bond can break easily due to the brittle structure of the material and should not be bent.
- A soft cleaning stone is recommended for the cleaning of tools with sintered bond. This procedure will remove the top layer of the bond and release unused diamond grit.

### Cutting discs with penetrated rim

In addition to proven cutting discs with surface layer for all grinding and cutting operations, we are the only manufacturer of cutting discs with penetrated rim covering the diameter range from 10 to 200 mm. This layer is - similar to a sintered bond - completely penetrated with diamond crystals. The discs allow thinnest cuts and extremely long service. They are constantly flexible and appropriate for technicians during cutting and lateral grinding of ceramics, porcelain, plaster and metals.

### Directions for use of rotary diamond discs with penetrated rim (Lab tools with electroplated bond, dispersed, Fig. 354).

Cutting discs with electroplated bond and penetrated rim provide unique properties. Other than discs with a single layer, this type achieves a very long life time by a specific multilayer design. Though they are thin and flexible, much more than the ones with sintered bond. They surpass all other types of discs when handled properly.

- The tools should be set to the working speed according to the graphical diagram.
- Cutting discs with penetrated rim are designed for separating. An oscillating motion should be applied while moving in radial direction.
- Although these discs are flexible, lateral grinding should be performed at very low forces only. High lateral forces may result in breakage, as the thickness of the blade usually is only 0,12 mm.
- Bending or turning when cutting into a slot must be avoided.
- A soft cleaning stone is recommended for the cleaning of tools with sintered bond. This procedure will remove the top layer of the bond and release unused diamond grit.

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
## Information and instructions for cleaning and sterilizing electroplated dental diamond burs with FG or right angle shank

- 1) *Dental Diamond burs with FG or angle shank made by fis can be reused in principle. Provided, however, that the tools will be processed before being used again. Burs that have already been used in grinding must be cleaned and sterilized. This similarly applies to brand new cutters. The latter must be sterilized before the first use.*
- 2) *This work should be solely carried out by employees who have acquired and demonstrated the necessary expertise and skills by training, education or practical activities.*
- 3) *Additional instructions are required in case of modifications or update in procedure or upon introduction of new products.*
- 4) *It should be ensured that the recommendations of Robert Koch Institute and Federal Institute for Pharmaceuticals and Medical Devices („Demands on hygiene for processing medical devices“) are followed. Special needs of dentistry should be taken into account.*
- 5) *The characteristics of rough surfaces of grinding tools impose special requirements on cleaning and sterilizing of those burs.*
- 6) *The grinding layers of dental diamond burs first need to be thoroughly rinsed under running water and cleaned with a stiff brush (not metal wire brush) in order to remove remains.*
- 7) *The burs should be subsequently cleaned in an ultrasonic bath at a temperature of 65° Celsius for about 1 to 3 minutes. Instructions of the device manufacturer should be observed. Grinding layers should not get in touch with the shanks of other tools while in the ultrasonic bath.*
- 8) *After cleaning in the ultrasonic bath the burs need to be inspected with a loupe (at least 6-fold). Step 5) and 6) should be repeated when indicated. In case of excessive wear, (partial) disbonding or sticky remains the burs should be disposed of.*
- 9) *All diamond burs with a functional layer should be then steam sterilized in a fractional vacuum process. The device should be designed according to EN 13060, grade B.*
- 10) *Taking into account the information and instructions of the manufacturer the following steps and settings should be selected:*
  - *fractional preliminary vacuum (4-fold)*
  - *sterilization temperature 134° Celsius*
  - *dwel time 5 minutes (full cycle)*
  - *drying time 10 minutes*
- 11) *After sterilization and drying the burs are to be packed in suitable packaging sterile.*
- 12) *Documented release and storage of sterilized and packaged dental diamond burs.*

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
Naturdiamantkristalle  
*Natural diamond crystals*

## Natural diamond - synthetic diamond

We manufacture our FG tools in varying grit sizes between 504 and 554 (acc. to ISO standards), with choice of one of two types of grit: natural and/or synthetic ("S") diamond (limited to ISO 524 and certain types). The natural diamond crystals typically have an irregular shape and sharp cutting edges which ensure a quick material removal. Layers containing the block-type, synthetic diamond crystals provide a professional preparation with less surface roughness. Those FG diamonds have a black coating to increase the optical contrast between tool and tooth.

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**Finzler, Schrock & Kimmel GmbH**